


Eliminate TFA and Improve Sensitivity of Peptide Analyses by LC–MS

TFA (trifluoroacetic acid) is a commonly used mobile phase additive for reversed-phase HPLC (RP-HPLC) separations of proteins and peptides. However, TFA interferes with and significantly reduces the LC–MS signal, lowering sensitivity. The ideal column for modern RP-LC–MS analysis should provide symmetrical peak shape without TFA in the mobile phase. The highly inert surface of Discovery BIO silica results in columns that give symmetrical and efficient peaks for peptides without TFA for maximum LC–MS sensitivity.

A major challenge facing biotechnology and proteomics researchers and others working with peptides or peptide maps is the need to detect and identify single peptides often at very low concentrations in extremely complex samples. Liquid chromatography coupled with mass spectroscopy (LC–MS) has become an invaluable tool to meet this challenge. Because of its high resolving power, RP-HPLC is the preferred separation mode for peptides. Traditionally, TFA is used in the mobile phases for RP-HPLC peptide separations. Ionic mobile phase additives such as TFA serve one or more of the following functions: pH control (buffering), complexation with oppositely charged ionic groups to enhance RP retention (ion pairing), or suppression of adverse ionic interactions between peptides and silanol groups on the silica. The latter function is necessary when using RP-HPLC phases with high silanol activity.

While TFA has little effect on UV detection, it has serious disadvantages for LC–MS detection. First, typical concentrations of TFA (0.1% v/v) have high surface tension and prevent efficient spray formation (nebulization). Second, TFA ions in the gas phase ion-pair with the peptide's basic groups suppressing their ionization and reducing sensitivity. A demonstration of TFA's adverse effect on LC–MS sensitivity is shown in Figure 1. Without TFA, the MS is able to detect much lower concentrations of these peptides. An added benefit is that at low TFA concentrations, resolution is improved because small differences in peptide retention are not masked. This is shown in the increased separation of peptides 1 and 2 in Figure 1 as the TFA concentration is decreased. At 0.1% TFA, they co-elute. Therefore,

Figure 1: Effect of chromatographic conditions on MS signals of peptides.



Column: Discovery BIO Wide Pore C18, 15 cm \times 2.1 mm, 3 μ m
Cat. No.: 567202-U
Mobile Phase: (A1) A: 25 mM formic acid in H₂O, B: 50:50 (25 mM formic acid in H₂O):(20 mM formic acid in CH₃CN)^a;
 (A2) A: 0.01% TFA, B: 0.01% TFA in 50:50 (CH₃CN:H₂O);
 (A3) A: 0.1% TFA, B: 0.1% TFA in 50:50 (CH₃CN:H₂O)
Flow-rate: 0.208 mL/min^b
Det.: +ES
Temp.: ambient
Inj.: 1 μ L or 3 μ L
Sample: RP Peptide Performance Standard, p/n RPS-P0010 (Alberta Peptide Institute)

a) molarity of formic acid adjusted to provide minimum baseline drift

b) linear velocity equal to 1 mL/min on 4.6 mm i.d. columns

Peptide 1: RGAGGLGLGK-amide

Peptide 2: ac-RGGGGLGLGK-amide

Peptide 3: ac-RGAGGLGLGK-amide

Peptide 4: ac-RGVGGLGLGK-amide

Peptide 5: ac-RGVVGLGLGK-amide

from the mobile phase standpoint, the best LC–MS method employs ionic additives other than TFA that are still volatile, can provide pH control, and do not strongly ion-pair with the peptides.

However, as previously mentioned, it is necessary to add TFA to the mobile phase when using RP-HPLC columns that exhibit high silanol activity. Without TFA, peptides, especially basic peptides, elute with low efficiency and tailing peaks and concurrently decreased sensitivity. Figure 2 shows a mixture of basic peptides on two columns that exhibit differing degrees of silanol activity. Without TFA, peptides on the popular brand column in Figure 2(b) show very asymmetrical, low efficiency peaks. However, because of its high degree of silanol deactivation and silica purity, the Discovery BIO wide-pore C18 column shown in Figure 2(a) provides efficient and symmetrical peaks. Therefore, from the column standpoint, for the best LC–MS method choose a column that has low surface silanol activity.

There are several TFA alternatives for LC–MS of proteins and peptides. For low pH operation, the most common are formic acid (HCOOH) and acetic acid (CH₃COOH). Ammonium acetate (CH₃COONH₄ or NH₄OAc) is commonly used for neutral pH

operation and ammonium bicarbonate (NH₄HCO₃) for basic pH. However, ionic additives not only control the pH, they can also influence the selectivity (peak order or spacing between peak apices). When choosing an ionic additive for LC–MS, it is important to consider its volatility, purity, pK_a, and degree to which it suppresses or expresses ionization of the analyte.


Of the myriad considerations for optimizing the LC–MS separation of proteins and peptides, two of the most important are to avoid TFA, and to use an RP-HPLC column that provides high efficiency and peak symmetry under those difficult conditions. Discovery BIO wide-pore columns and capillaries provide sensitive, efficient, stable and reproducible LC–MS analyses of proteins and peptides without TFA in the mobile phase.

Figure 2: Toward basic peptides without TFA.

(a) Discovery BIO Wide Pore C18, 300 Å, 15 cm × 2.1 mm i.d., 5 µm, 0.208 mL/min

Gradient


Min	°C	%D
0	80	20
20	40	60



(b) Competitor C18, 300 Å, 15 cm × 2.0 mm i.d., 5 µm, 0.189 mL/min

Gradient

Min	°C	%D
0	80	20
20	40	60



Columns: C18, 300 Å, 15 cm × 2.1 mm or 2.0 mm i.d., 5 µm

Mobile Phase: A: 25 mM formic acid in water
B: 50:50 (25 mM formic acid in water):(20 mM formic acid in CH₃CN)

Flow Rate: 0.208 (or 0.189) mL/min

Det.: 215 nm

Temp.: 35 °C

Inj.: 0.5 µL (~0.25 µg ea peptide)


Sample: RP Peptide Ionic Interactions Standard,
p/n RPS-I0020
(Alberta Peptide Institute)

Peptide 1: ac-GGGLGGAGGLK-amide

Peptide 2: ac-KYGLGGAGGLK-amide

Peptide 3: ac-GGALKALKGLK-amide

Peptide 4: ac-KYALKALKGLK-amide



SUPELCO

Sigma-Aldrich Company Ltd/Supelco UK,
Fancy Road, Poole, Dorset BH12 4QH, UK,
tel. +44 1202 716789,
fax +44 1202 715460,
e-mail: ukorders@eurnotes.sial.com
website: www.sigma-aldrich.com/supelco